

How do I know which face mask to use?

Order	Droplet precautions (not COVID)	"Universal masking" for pandemic		
		Healthcare worker		Patient and visitor
		Direct patient care	Not direct patient care	
#1	ASTM Level 1	ASTM Level 1	ASTM Level 1	ASTM Level 1
#2	ASTM Level 1	ASTM Level 1	Type B	Type B
#3	Type A	Type A	Type B	Type B
#4	Type A	Type A	Type B	None
#5	Type B	Type B	Type B	None
#6	Cloth masks	Cloth Masks	Cloth Masks	None
			Current products in use	

Notes

- Type A and Type B masks have all been assessed as acceptable for use for the above indications
- N95 and N95-equivalents are to be used in accordance with latest PPE guidance.
- ASTM Level 2 & 3 face masks are to be used only with procedures with high risk of fluid splash in settings such as operating rooms and IR/cath areas
- For available product pictures and INFOR details refer to the COVID-19 Intranet Site - <https://www.ynhhs.org/patient-care/covid-19/for-employees/ppe.aspx>

How do I know which protective clothing to use?


Order	Droplet precautions (not COVID)
#1	Isolation Gowns
#2	<p><u>Alternative Protective Clothing</u></p> <ul style="list-style-type: none"> • Lab Coats and Frocks <ul style="list-style-type: none"> • Coveralls • Aprons • Other

Current product in use

Notes

- Products in the current status or higher are acceptable for use
- Surgical gowns are used only for sterile procedures (e.g. operating room)
- For available product pictures and INFOR details refer to the COVID-19 Intranet Site - PPE link

How do I know which fit tested respirator (FFR) to use?

	I normally use: 	3M 1860 Small	3M 1860	3M 1870	Moldex
	#1	If you have been fit tested for an elastomeric respirator, this is your primary fit respirator			
Order of use	#2	3M 1860S	3M 1860	3M 1870	Moldex XS
	#3	3M 8110S	3M 8210, 3M 8210 plus	3M 9205 3M 9210	Moldex XS
	#4	PAPRS are used if fit tested respirators are not available or you have not been able to be fit tested to any of the respirators above			
	#5	Unfit tested FFRs			
	#6	Facemask (will be used only if there are no remaining FFRs – this is unlikely to occur)			

Current product in use

- | | |
|-------|---|
| Notes | <ul style="list-style-type: none"> • Products in the current status or higher are acceptable for use • #1 provide similar fit as your usual respirator but may be industrial and require additional PPE if used during while in the sterile field during a procedure • Fit tested respirators should only be used in accordance with the latest PPE guidance. • For available product pictures and INFOR details refer to the COVID-19 Intranet Site - PPE link |
|-------|---|