Ca	ategory	Product Image	Product #	INFOR#	Product Name	Manufacturer Name	Manufacturer Product #	NOTES	Use with a face mask in OR/sterile procedure (Yes/No)													
			17-986-260	220622	EVA PAPR SYSTEM W/LOOSE FITTING FACEPIECE LG	Bullard	EVA20LF2L	System comes with 1 large hood and 1 breathing tube	Not applicable		3		E	ntire u	nit for use inlo	cudes H	ood, Blo	wer unit	t, Hose,E	Battery a	nd Cha	rger
			18-998-002	220627	LF2L SERIES LOOSE FIT HOOD LARGE	Bullard	20LF2L	More of a narrow fit	Yes													
			18-998-003	220628	LF2M SERIES LOOSE FIT HOOD MEDIUM	Bullard	20LF2M	More of a narrow fit	Yes													
			19-021-233C	220620	LFL SERIES LOOSE FIT LARGE HOOD	Bullard	20LFL	Loose Fit Large hood	Yes													
	P A P R		19-098-043	221234	LF2L SERIES WITH HEPA FILTER BUILT IN	Bullard	20LF2LHE	On backorder	Yes													
	.,	00	17-986-282	220621	FILTER CARTRIDGE EVA PAPR HEPA FILTER RESPIRATOR	Bullard	PAPRFC3	3 Week usage before replacement	Not applicable													
			17-986-277	220624	BATTERY PACK EVA PAPR BATTERY 7-8 HR	Bullard	EVABAT1	10-12 Hour life at full battery. 4 hour charge to full battery.	Not applicable													
	ľ	0.	19-021-238	220619	BREATHING TUBE REPLACEMENT	Bullard	PA20LFBT	Replace when integrity of tube is damaged.	Not applicable													
		X	NC1716047	220674	BACKPACK SUSPENDERS	Bullard	EVABKPK2	Assist with distributing weight and holding up belt	Not applicable													
		Bullers	NC1810720	220675	DUFFLE BAG STORAGE	Bullard	PAPRBAG	Fits 2 PAPR units	Not applicable													
			17-986-281	220625	GANG CHARGER BATTERY PACK	Bullard	EVAGC	Charge up to 6 batteries at one time	Not applicable													
	p																					
							•													'	'	

	Product Image	Product #	INFOR#	Product Name	Manufacturer Name	Manufacturer Product #	NOTES	Use with a face mask in OR/sterile procedure (Yes/No)							
A P															
R															



Entire Unit includes all seven Items listed

Category	Product Image	Product #	INFOR#	Product Name	Manufacturer Name	Manufacturer Product #	NOTES	Use with a face mask in OR/sterile procedure (Yes/No)
				PAPR unit	ILC Dover (Sentinel)			Not applicable
	The state of the s			Head Cover	ILC Dover (Sentinel)	S-4001		Yes
				CBRNE Hood	ILC Dover (Sentinel)			No
				CBRNE Cartridge		P.N. S-2016		Not applicable
				BREATHING TUBE REPLACEMENT	ILC Dover (Sentinel)			Not applicable
	Name of the first of the second of the secon			D Cell battery adapter kit	ILC Dover (Sentinel)	S-4013		Not applicable